

## NEEDLE KIT IDENTIFICATION PACKAGE



- Unique and permanent pet identification
- Bio-compatible glass transponder, according to ISO standard 11784/5
- Easy and quick microchip implantation
- Reliable proof of ownership
- Single use needles, treated by EtO gas
- Implanter made in resistant and light plastic material
- Ideal product for city or government tenders
- ISO-compliant microchips can be read nationally and internationally

### Product description and benefits:

#### Transponder:

Datamars new Needle kit gives pets a **silent but effective voice**. Unlike a collar with tag or tattoo, a microchip cannot be falsified, removed, torn or chewed off. Each glass transponder is as **small as a grain of rice** (13.3 x 2.12 mm), and contains a permanent, unalterable ID **code that is unique to that animal**.

The microchip comes pre-loaded in a disposable and single-use needle, **ready for use in combination with the provided pistol-type applicator**. The applicator has a male luer-lock connection fitting to be fit together with the needle, which screws into it by twisting it clockwise. The ergonomics of the red plastic implanter, extremely light (only 42 g), makes the implantation of the microchips easy and quick. Each blister includes individually barcode labels to simplify documentation and record keeping and has coloured indication to confirm that it has been **treated by EtO gas**.


The Needle kit product, thanks to its combination of quality and competitive cost, is **best suited for city or government tenders** that have a need for larger microchip quantities. Once implantation is complete, the detachable **needle can be removed for easy disposal**.

## RECOMMENDED IMPLANTATION SITES

**Canine and Feline:** The transponder is implanted subcutaneously on the left side of the neck, behind the ear, lateral of the fourth to the fifth cervical vertebra.


**Equine:** The microchip is implanted within the nuchal ligament in its middle third or at the halfway point between the ears and the withers. This is the recommended implantation site in all countries except Australia where the microchip is implanted in the musculature of the left neck or the anterior injection triangle.

**For more information on recommended microchip implantation sites, visit the World Small Animal Veterinary Association at [www.wsava.org](http://www.wsava.org).**



### Reader:

The reader communicates with RFID transponders via low power, low frequency radio waves. The passive transponder transmits its preprogrammed unalterable code, identifying the animal at high speed. The transponder can be read through any non-conductive material.



### KEYS TO EFFECTIVE SCANNING:

1. Hold the scanner close to or touching the pet
2. Scan slowly and patiently
3. Rock the scanner back and forth slightly while scanning because implanted chips are in various orientations
4. Begin and concentrate scanning at the cranio-dorsal aspect of the pet
5. Scan in a horizontal then a vertical "S" pattern down the pet as shown in the diagram
6. Consider scanning each pet more than once

### CORRECT IMPLANTATION PROCEDURE:

1. Scan the pet to check if it carries a transponder already
2. Verify function and correct number of the transponder prior to implantation
3. Implant the transponder according to recommended implantation sites. After implantation, scan the microchip again to verify function and don't forget to register it in a pet recovery database
4. Twist the needle assembly anticlock-wise to remove it from the body of the syringe and dispose it according to local regulations for veterinary waste



### TECHNICAL DATA:

Size (L x Ø)	13.3 mm ± 0.4 x 2.12 mm ± 0.05
Weight	≈0.114 g
Operating temperature	-25°C to +70°C
Storage temperature	-40°C to +90°C
Passive, without battery	
ID code	15 digits conforming to ISO Standard 11784, FDX-B type
Working Frequency	134.2 kHz
Memory capacity	Total memory 512 bits, 64 bits required for ISO Standard 11784
Reading distance	Up to 30 cm
Housing	Bio-compatible glass
Packaging	Box of 100 transponders each

Subject to change without notice



ICAR MEMBER

DATAMARS' innovative RFID-based solutions are the result of more than 25 years of industry experience that have enabled the company to establish itself as the global market leader in the companion animal identification market. Having driven the development of ISO standards for companion animal identification, DATAMARS works with pharmaceutical companies, shelters, kennel clubs and charitable organizations to deploy its unique pet identification and reunification solution to improve the lives of companion animals around the world.

## DATAMARS

### Americas Headquarters:

250 West Cummings Park  
Woburn, MA 01801, USA  
Phone: +1 781 281 2216  
Fax: +1 781 300 7330

### Corporate Headquarters:

Via ai Prati  
6930 Bedano-Lugano  
Switzerland  
Phone: +41 91 935 73 80  
Fax: +41 91 945 03 30

### Asia Headquarters:

Northern Region Industrial Estate  
179/1 Moo 4, T. Ban Klang, A. Muang  
Lamphun, Thailand 51000  
Phone: +66 53 582 021  
Fax: +66 53 582 680

[animal-id@datamars.com](mailto:animal-id@datamars.com)

[www.datamars.com](http://www.datamars.com)