


S-ID Software

Datamars S-ID software is designed to work specifically with Datamars companion animal and livestock readers to help simplify data transfer and reader management associated with electronic animal identification.

S-ID sports a simple, intuitive interface with a wizard to help streamline the setup process. S-ID will automatically detect which Datamars reader is attached and will then display the appropriate functions and options. S-ID will help with data transfer to computer or application, and will also be the interface for downloadable updates and upgrades.



Benefits

- Simple, intuitive interface makes S-ID easy to use out-of-the-box
- Simplifies reader update and upgrade process
- Provides easy to use data management tools and helps with data organization
- Connects to the PetLink database and Petmaxx search engine

Features

- Automatically checks for software and reader updates
- Data management functions, including browsing, downloading, editing, export of .csv and .xml, and data deletion
- Backup and restore function
- Quick connection with Datamars pet databases and search engines worldwide
- Can be used to send remote commands to the reader (ex.: remote start of microchip reading)