

COMPACT MAX

Portable Reader for Animal Identification

COMPACT MAX reads FDX-B (ISO) and FDX-A/Fecava microchips, featuring best in class reading distances in a pocket size dimension.

COMPACT MAX has an ergonomic design and an intuitive icon-based interface, making the reader perfect for animal identification, from puppies and kittens to seniors.

COMPACT MAX can read up to 4'000 microchip codes and can easily transfer data to your PC via USB.



Highlights

- Best-in-class reading performances, able to read FDX-B and FDX-A/Fecava microchips
- Pocket size and ergonomic design
- Bright color display
- Up to 4'000 microchip codes stored in the internal memory
- Easy data transfer to a PC

Dimensions	16 x 80 x 2 cm	Interfaces	USB (Mass storage, VCP and Keyboard wedge)
Weight	170 g	Transmission frequency	134.2 KHz
Housing material	ABS	Memory	Up to 4'000 codes
Operating temperature	From -5°C to +55°C, 95% RH, non condensing	Typical reading distance	FDX-B: up to 16 cm FDX-A/FECAVA: up to 8 cm
Charging temperature	From +5°C to +45°C, 95% RH, non condensing	Certifications	CE, FCC
Storage temperature	From -20°C to +65°C, 95% RH, non condensing	Standards	ISO 11784/5
Battery	Li-ion polymer 3.7V (1250 mAh)		

Subject to change without notice - April 2023

Datamars SA
Corporate Headquarters
Via Industria 16 6814 Lamone, Switzerland
Phone +41 91 935 73 80 | Fax +41 91 945 03 30
animal-id@datamars.com | pet.datamars.com

DATAMARS
Animal iD