

COMPACT MAX+ Universal Scanner



COMPACT MAX+ is a low frequency microchip scanner able to read FDX-B in accordance with the ISO Standard 11784/5, FDX-A, Trovan and Avid. COMPACT MAX+ features the best in class reading distances in a pocket size dimension and can store up to 4,000 microchip numbers in its internal memory. This scanner has an integrated low energy Bluetooth module, which allows the reader to have all the advantages of a wireless connection while providing long-lasting battery life.

Highlights

- Integrated low energy Bluetooth module
- Long-lasting rechargeable Li-ion battery
- Bright color display
- Universal, reads ISO FDX-B, FDX-A, Trovan and Avid
- Dedicated app to manage stored data and advanced settings

Technical Specifications

Dimensions	3.14" x 6.2" x 0.7" in	Interfaces	Bluetooth 5.2 (BLE), USB Type-C (Mass storage, VCP and Keyboard Wedge)
Weight	6.07 oz	Transmission frequency	134.2 KHz
Housing material	ABS	Memory	Up to 4,000 microchip numbers
Operating temperature	From 23°F to 131°F, 95% RH, non condensing	Typical reading distance	FDX-B: up to 7.09" in / FDX-A / Fecava: up to 3.54" in AVID ENCRYPTED: up to 3.94" in TROVAN: up to 2.17" in
Charging temperature	From 41°F to 113°F, 95% RH, non condensing	Certifications	CE, FCC, UKCA
Storage temperature	From -4°F to 149°F, 95% RH, non condensing	Standards	ISO 11784/5
Battery	Li-ion polymer 3.7V (1800 mAh)		

PetLink Mailing Address:

Datamars, Inc. 345 West Cummings Park, Woburn, MA 01801

Contact us at

1-877-PETLINK (738-5465)

or email us at

petlink@petlink.net

www.petlink.net

ussales@datamars.com



© 2025 Datamars, Inc. All rights reserved. PetLink is a registered trademark of Datamars, Inc.